

# PUFFY EYE TREATMENT PROTOCOL

**Step 1:**

Ultrasound Cleansing with HydroClear Serum  
Step duration: 20-25 minutes

**Step 1:**

Vacuum Therapy with DazzleSkin Pure Plant-Based Nourishing Massage Oil  
to improved blood and lymphatic circulation  
Step duration: 10-15 minutes

**Step 3:**

Electroporation with WrinkleGone Peptide Rejuvenation Serum  
Step duration: 10-15 minutes

**Step 4:**

Cooling therapy with FirmingGenius Peptide Conductive gel  
Step duration: 10 minutes

**Step 5:**

Apply moisturizer and sunscreen SPF 40.

